

**PATENT**

**APPLICATION FOR UNITED STATES LETTERS PATENT**

**for**

**CLARIFICATION OF PROTEIN PRECIPITATE SUSPENSIONS  
USING ANIONIC POLYMERIC FLOCCULANTS**

**by**

**MARK L. CHIVERS**

Express Mail Mailing Label Number:  
Date of Deposit:

EL291388044US  
July 16, 1999

0933931 032704  
T0220 T06550

## BACKGROUND OF THE INVENTION

This is a conversion of provisional application serial number 60/093,555 filed July 21, 1998.

The present invention relates generally to the purification of proteins produced by recombinant DNA technologies. More particularly, it concerns the clarification of protein suspensions containing soluble and insoluble components via flocculation with anionic polymers.

The use of flocculating agents has been described in several industrial settings, including the biotechnology industry. For example, U.S. Patent No. 5,047,511 describes the use of cationic flocculating agents in the recovery of recombinant somatotropin protein from a protein solution containing soluble high molecular weight contaminating proteins. This involved the selective precipitation of contaminating high molecular weight proteins by adding a cationic flocculant containing quaternary ammonium groups and then recovering the somatotropin from the solution.

Somatotropins, also known as growth hormones, are polypeptides produced and secreted by cells of the pituitary gland. These proteins are known to be effective in promoting pre-adult skeletal growth and meat production of beef cattle and swine, and can be produced reliably and inexpensively in large quantities by recombinant DNA technology. In addition, they are known to affect a variety of metabolic processes including the stimulation of lactation, improvements of the efficiency of converting feed to meat or milk, lipid-mobilizing effects, and others.

Recombinant DNA technology provides a means for the large scale production of heterologous proteins of interest in bacterial host cells. In the case of somatotropin, a growth hormone, the protein is sequestered in inclusion bodies within the cytoplasm of the host cells. The inclusion bodies can be recovered from the host cell culture by disrupting the cell so as to release the inclusion bodies, and thereafter collecting the inclusion bodies as a solid pellet by differential centrifugation. The inclusion bodies are solubilized in an aqueous solution of a suitable chaotropic agent such as urea or guanidine hydrochloride at an alkaline pH and subsequently naturized by contact with a mild oxidizing agent to form intramolecular disulfide

bonds and to refold the protein into its biologically active, native conformation. Methods for the solubilization and naturation of somatotropin protein produced by *E.coli* bacteria are described in U.S. Patent No. 4,511,502 and U.S. Patent No. 4,652,630, each of which is incorporated herein by reference.

5           The somatotropin refold solution obtained from the naturation step (as described for example in U.S. Patent Nos. 4,511,502 and 4,652,630) comprises an aqueous solution of somatotropin monomers, dimers and higher oligomers, along with residue and other debris from the host cells. Of these, the somatotropin monomer is the desired biologically active agent. U.S. Patent No. 5,182,369, the disclosure of which is incorporated herein by reference, describes the  
10       selective precipitation of somatotropin dimer and higher oligomers together with residual host cell proteins and other contaminating substances from a pH-adjusted somatotropin refold solution, leaving the desired somatotropin monomer as the primary soluble constituent of the suspension.

          Once the somatotropin oligomers and other contaminants have been selectively  
15       precipitated using this approach, it is necessary to remove the precipitated proteins and other insoluble contaminants from the suspension in order to obtain somatotropin monomers of the desired purity. Such liquid/solid separations as those required for this purification step are employed in most industrial biotechnological processes and are frequently accomplished via centrifugation and/or filtration procedures.

20       Flocculating agents can be employed to improve liquid/solid separations by aggregating the solids that are present in a protein suspension, thereby increasing the particle size of the solids (for review, see Halverson and Panzer, 1980). An increase in particle size is particularly beneficial in centrifugation and sedimentation applications where the particle sedimentation velocity is proportional to the square of the particle radius. The increased sedimentation velocity  
25       that results from larger particle sizes can improve productivity in any type of liquid/solid separation where particle sedimentation velocity is a factor.

## SUMMARY OF THE INVENTION

This invention broadly concerns the separation of soluble proteins from insoluble contaminants via flocculation. More particularly, it relates to the use of anionic polymeric flocculants in the separation and recovery of somatotropin proteins.

5 Therefore, in accordance with one aspect of the present invention, there is provided a method for separating an aqueous protein suspension of soluble somatotropin monomer and insoluble contaminants by adding to the suspension an anionic polymer in an amount and under conditions effective to cause the flocculation of the insoluble contaminants. The flocculated, insoluble, material can be easily separated from the soluble somatotropin monomer to recover a  
10 clarified supernatant of soluble somatotropin monomer.

In accordance with another aspect of the invention, there is provided a method for the recovery of somatotropin monomer which comprises

obtaining a mixture of somatotropin proteins comprising somatotropin monomer and somatotropin oligomer in aqueous solution at a pH greater than about 7;

15 producing a protein suspension containing soluble somatotropin monomer by precipitating a major portion of the somatotropin oligomer from the solution while maintaining a major portion of somatotropin monomer in solution by reducing the pH of the solution to less than about 6.5;

20 adding to the protein suspension containing soluble somatotropin monomer an anionic polymer in an amount and under conditions effective to cause the flocculation of the precipitated proteins; and

separating the flocculated material from the solution of somatotropin monomer.

25 In accordance with another aspect of the invention, there is provided an aqueous protein suspension comprising somatotropin monomers, somatotropin oligomers (oligomer as used herein refers to dimers as well as other multimeric forms of the protein), and an anionic polymer.

Suitable anionic polymers used in accordance with the method of this invention include but are not limited to polyacrylamides, particularly those having charge densities in the range of about 1 to about 20%, and polysaccharides such as starch.

5 This invention provides an improved means for separating somatotropin monomer from the insoluble contaminants that are selectively precipitated from the monomer during recovery of the recombinant protein. Flocculation of the insoluble contaminants increases their particle sizes and therefore their sedimentation velocities, providing improved productivity in liquid/solid separations.

### 10 BRIEF DESCRIPTION OF THE DRAWINGS

The following drawings form part of the present specification and are included to further demonstrate certain aspects of the present invention. The invention may be better understood by reference to one or more of these drawings in combination with the detailed description of specific embodiments presented herein.

15 FIGURE 1 illustrates the effect of polymer charge density on the flocculation performance.

FIGURE 2 illustrates the effect of pH on the flocculation performance from about pH 4.0 to about pH 6.0.

FIGURE 3 illustrates the effect of pH on the flocculation performance from about pH 4.0 to about pH 4.8.

### 20 DESCRIPTION OF ILLUSTRATIVE EMBODIMENTS

This invention provides anionic polymeric flocculants for improving the efficiency of liquid/solid separations by increasing the particle size of insoluble contaminants in a protein suspension.

25 "Protein suspension," as used herein, represents an aqueous suspension, preferably at an acidic pH (i.e., less than about 7), which comprises soluble and insoluble components.

Typically, but not necessarily, the soluble component includes the desired protein product (e.g. soluble somatotropin monomer) whereas the insoluble components include the undesired contaminants (e.g. somatotropin oligomers and bacterial cell debris).

“Somatotropin” or “ST” refers to any polypeptide that has biological activity and/or chemical structure similar to that of a somatotropin produced in the pituitary gland of an animal. Such somatotropins include natural somatotropins produced by the pituitary somatotrophic cells or, alternatively, somatotropins produced by recombinant DNA technology in which a somatotropin or a variant derived therefrom is expressed by genetically transformed prokaryotic or eukaryotic cells, e.g. bacteria (such as *E. coli*), yeast or algae.

“Refold solution” refers to the stock solution obtained as a result of the folding and oxidation in the somatotropin naturation step, as described for example in U.S. Patent Nos. 4,511,502 and 4,652,630. Subsequent adjustment of the pH of the refold solution to below about 6.5, preferably to between about 4.0 and about 6.0, selectively precipitates the insoluble contaminants while leaving the soluble somatotropin monomer substantially unaffected.

“ST protein suspension” refers to the aqueous suspension produced essentially in accordance with U.S. Patent No. 5,182,369, such that the suspension has a pH between about 4.0 and about 6.5 and comprises primarily soluble somatotropin monomer and insoluble contaminants. Prior to the pH precipitation, the solution generally contains from about 1 to about 30 grams/liter total protein comprising somatotropin, somatotropin aggregates and *E. coli* proteins. Approximately 40-60% of the total protein is typically precipitated in the pH adjustment step, resulting in an oxidized somatotropin monomer purity greater than about 90% in the liquid fraction.

“Flocculation” refers to the aggregation of insoluble particles caused by the addition of a suitable flocculating agent to a protein suspension. By increasing the particle size of the insoluble components present in the suspension, the efficiency of solid/liquid separations, such as by filtration or centrifugation, is improved.

“Clarification” refers to the removal of insoluble contaminants from a protein suspension, thereby improving the clarity/purity of the suspension, for example as measured by a decrease in the spectrophotometric absorbance of the suspension at 700 nm.

5 “Separation” refers to the removal of the aqueous phase, containing the soluble somatotropin monomers, from the insoluble particles following flocculation of the ST protein suspension. This removal is accomplished by any means compatible with the present invention including the common industrial methods such as sedimentation and filtration. This separation results in the recovery of a solution of somatotropin monomers which is essentially free of insoluble contaminants.

10 “Sedimentation” refers to the settling of the flocculated precipitated contaminants, either by centrifugation or by gravity.

The following art, to the extent that it provides exemplary procedural or other details supplementary to those set forth herein, is specifically incorporated herein by reference:

15 Garcia, F.A.P., “Protein Precipitation” in: Recovery of Biological Materials, John Wiley & Sons, 1993.

Gates et al., “Selecting Agitator Systems to Suspend Solids in Liquids”, Chemical Engineering, May 24, 1976.

20 Halverson, F., Panzer, H.P., “Flocculating Agents” in: Encyclopedia of Chemical Technology, Volume 10, Third Edition, John Wiley & Sons, 1980, pp 489-523.

Muhle, K. and Domasch, K. “Stability of Particle Aggregates in Flocculation with Polymers”, *Chemical Engineering Progress* 29 (1991).

Wurzburg O.B., “Modified Starches: Properties and Uses”, CRC Press, Inc., pp. 4-98.

25 In one aspect of this invention, there are provided methods for improved liquid/solid separations of an aqueous ST protein suspension containing soluble proteinaceous product and insoluble contaminants. This method involves the clarification of a protein suspension by

10  
5  
floculation of the precipitated contaminants using anionic polymers. A suitable flocculant is dispersed in a protein suspension under conditions effective to cause the aggregation of the insoluble proteins and other contaminants present in the suspension while leaving the soluble product substantially unaffected. A number of advantages are provided by the flocculation procedure described herein. First, separation efficiency is improved due to the fact that the sedimentation velocity of solids increases with increasing particle size. This can effectively minimize or even eliminate the need for centrifugation in a solid/liquid separation, thereby reducing the capital costs associated with the process. Also, the method can improve filtration processes by reducing the pressure drop that is often observed across filters. This change in hydrodynamic properties can allow for higher flow rates and/or the use of less diatomaceous earth (i.e., filter aid) in filtration processes.

15  
20  
25  
30  
35  
40  
45  
50  
55  
60  
65  
70  
75  
80  
85  
90  
95  
100  
105  
110  
115  
120  
125  
130  
135  
140  
145  
150  
155  
160  
165  
170  
175  
180  
185  
190  
195  
200  
205  
210  
215  
220  
225  
230  
235  
240  
245  
250  
255  
260  
265  
270  
275  
280  
285  
290  
295  
300  
305  
310  
315  
320  
325  
330  
335  
340  
345  
350  
355  
360  
365  
370  
375  
380  
385  
390  
395  
400  
405  
410  
415  
420  
425  
430  
435  
440  
445  
450  
455  
460  
465  
470  
475  
480  
485  
490  
495  
500  
505  
510  
515  
520  
525  
530  
535  
540  
545  
550  
555  
560  
565  
570  
575  
580  
585  
590  
595  
600  
605  
610  
615  
620  
625  
630  
635  
640  
645  
650  
655  
660  
665  
670  
675  
680  
685  
690  
695  
700  
705  
710  
715  
720  
725  
730  
735  
740  
745  
750  
755  
760  
765  
770  
775  
780  
785  
790  
795  
800  
805  
810  
815  
820  
825  
830  
835  
840  
845  
850  
855  
860  
865  
870  
875  
880  
885  
890  
895  
900  
905  
910  
915  
920  
925  
930  
935  
940  
945  
950  
955  
960  
965  
970  
975  
980  
985  
990  
995  
1000  
1005  
1010  
1015  
1020  
1025  
1030  
1035  
1040  
1045  
1050  
1055  
1060  
1065  
1070  
1075  
1080  
1085  
1090  
1095  
1100  
1105  
1110  
1115  
1120  
1125  
1130  
1135  
1140  
1145  
1150  
1155  
1160  
1165  
1170  
1175  
1180  
1185  
1190  
1195  
1200  
1205  
1210  
1215  
1220  
1225  
1230  
1235  
1240  
1245  
1250  
1255  
1260  
1265  
1270  
1275  
1280  
1285  
1290  
1295  
1300  
1305  
1310  
1315  
1320  
1325  
1330  
1335  
1340  
1345  
1350  
1355  
1360  
1365  
1370  
1375  
1380  
1385  
1390  
1395  
1400  
1405  
1410  
1415  
1420  
1425  
1430  
1435  
1440  
1445  
1450  
1455  
1460  
1465  
1470  
1475  
1480  
1485  
1490  
1495  
1500  
1505  
1510  
1515  
1520  
1525  
1530  
1535  
1540  
1545  
1550  
1555  
1560  
1565  
1570  
1575  
1580  
1585  
1590  
1595  
1600  
1605  
1610  
1615  
1620  
1625  
1630  
1635  
1640  
1645  
1650  
1655  
1660  
1665  
1670  
1675  
1680  
1685  
1690  
1695  
1700  
1705  
1710  
1715  
1720  
1725  
1730  
1735  
1740  
1745  
1750  
1755  
1760  
1765  
1770  
1775  
1780  
1785  
1790  
1795  
1800  
1805  
1810  
1815  
1820  
1825  
1830  
1835  
1840  
1845  
1850  
1855  
1860  
1865  
1870  
1875  
1880  
1885  
1890  
1895  
1900  
1905  
1910  
1915  
1920  
1925  
1930  
1935  
1940  
1945  
1950  
1955  
1960  
1965  
1970  
1975  
1980  
1985  
1990  
1995  
2000  
2005  
2010  
2015  
2020  
2025  
2030  
2035  
2040  
2045  
2050  
2055  
2060  
2065  
2070  
2075  
2080  
2085  
2090  
2095  
2100  
2105  
2110  
2115  
2120  
2125  
2130  
2135  
2140  
2145  
2150  
2155  
2160  
2165  
2170  
2175  
2180  
2185  
2190  
2195  
2200  
2205  
2210  
2215  
2220  
2225  
2230  
2235  
2240  
2245  
2250  
2255  
2260  
2265  
2270  
2275  
2280  
2285  
2290  
2295  
2300  
2305  
2310  
2315  
2320  
2325  
2330  
2335  
2340  
2345  
2350  
2355  
2360  
2365  
2370  
2375  
2380  
2385  
2390  
2395  
2400  
2405  
2410  
2415  
2420  
2425  
2430  
2435  
2440  
2445  
2450  
2455  
2460  
2465  
2470  
2475  
2480  
2485  
2490  
2495  
2500  
2505  
2510  
2515  
2520  
2525  
2530  
2535  
2540  
2545  
2550  
2555  
2560  
2565  
2570  
2575  
2580  
2585  
2590  
2595  
2600  
2605  
2610  
2615  
2620  
2625  
2630  
2635  
2640  
2645  
2650  
2655  
2660  
2665  
2670  
2675  
2680  
2685  
2690  
2695  
2700  
2705  
2710  
2715  
2720  
2725  
2730  
2735  
2740  
2745  
2750  
2755  
2760  
2765  
2770  
2775  
2780  
2785  
2790  
2795  
2800  
2805  
2810  
2815  
2820  
2825  
2830  
2835  
2840  
2845  
2850  
2855  
2860  
2865  
2870  
2875  
2880  
2885  
2890  
2895  
2900  
2905  
2910  
2915  
2920  
2925  
2930  
2935  
2940  
2945  
2950  
2955  
2960  
2965  
2970  
2975  
2980  
2985  
2990  
2995  
3000  
3005  
3010  
3015  
3020  
3025  
3030  
3035  
3040  
3045  
3050  
3055  
3060  
3065  
3070  
3075  
3080  
3085  
3090  
3095  
3100  
3105  
3110  
3115  
3120  
3125  
3130  
3135  
3140  
3145  
3150  
3155  
3160  
3165  
3170  
3175  
3180  
3185  
3190  
3195  
3200  
3205  
3210  
3215  
3220  
3225  
3230  
3235  
3240  
3245  
3250  
3255  
3260  
3265  
3270  
3275  
3280  
3285  
3290  
3295  
3300  
3305  
3310  
3315  
3320  
3325  
3330  
3335  
3340  
3345  
3350  
3355  
3360  
3365  
3370  
3375  
3380  
3385  
3390  
3395  
3400  
3405  
3410  
3415  
3420  
3425  
3430  
3435  
3440  
3445  
3450  
3455  
3460  
3465  
3470  
3475  
3480  
3485  
3490  
3495  
3500  
3505  
3510  
3515  
3520  
3525  
3530  
3535  
3540  
3545  
3550  
3555  
3560  
3565  
3570  
3575  
3580  
3585  
3590  
3595  
3600  
3605  
3610  
3615  
3620  
3625  
3630  
3635  
3640  
3645  
3650  
3655  
3660  
3665  
3670  
3675  
3680  
3685  
3690  
3695  
3700  
3705  
3710  
3715  
3720  
3725  
3730  
3735  
3740  
3745  
3750  
3755  
3760  
3765  
3770  
3775  
3780  
3785  
3790  
3795  
3800  
3805  
3810  
3815  
3820  
3825  
3830  
3835  
3840  
3845  
3850  
3855  
3860  
3865  
3870  
3875  
3880  
3885  
3890  
3895  
3900  
3905  
3910  
3915  
3920  
3925  
3930  
3935  
3940  
3945  
3950  
3955  
3960  
3965  
3970  
3975  
3980  
3985  
3990  
3995  
4000  
4005  
4010  
4015  
4020  
4025  
4030  
4035  
4040  
4045  
4050  
4055  
4060  
4065  
4070  
4075  
4080  
4085  
4090  
4095  
4100  
4105  
4110  
4115  
4120  
4125  
4130  
4135  
4140  
4145  
4150  
4155  
4160  
4165  
4170  
4175  
4180  
4185  
4190  
4195  
4200  
4205  
4210  
4215  
4220  
4225  
4230  
4235  
4240  
4245  
4250  
4255  
4260  
4265  
4270  
4275  
4280  
4285  
4290  
4295  
4300  
4305  
4310  
4315  
4320  
4325  
4330  
4335  
4340  
4345  
4350  
4355  
4360  
4365  
4370  
4375  
4380  
4385  
4390  
4395  
4400  
4405  
4410  
4415  
4420  
4425  
4430  
4435  
4440  
4445  
4450  
4455  
4460  
4465  
4470  
4475  
4480  
4485  
4490  
4495  
4500  
4505  
4510  
4515  
4520  
4525  
4530  
4535  
4540  
4545  
4550  
4555  
4560  
4565  
4570  
4575  
4580  
4585  
4590  
4595  
4600  
4605  
4610  
4615  
4620  
4625  
4630  
4635  
4640  
4645  
4650  
4655  
4660  
4665  
4670  
4675  
4680  
4685  
4690  
4695  
4700  
4705  
4710  
4715  
4720  
4725  
4730  
4735  
4740  
4745  
4750  
4755  
4760  
4765  
4770  
4775  
4780  
4785  
4790  
4795  
4800  
4805  
4810  
4815  
4820  
4825  
4830  
4835  
4840  
4845  
4850  
4855  
4860  
4865  
4870  
4875  
4880  
4885  
4890  
4895  
4900  
4905  
4910  
4915  
4920  
4925  
4930  
4935  
4940  
4945  
4950  
4955  
4960  
4965  
4970  
4975  
4980  
4985  
4990  
4995  
5000  
5005  
5010  
5015  
5020  
5025  
5030  
5035  
5040  
5045  
5050  
5055  
5060  
5065  
5070  
5075  
5080  
5085  
5090  
5095  
5100  
5105  
5110  
5115  
5120  
5125  
5130  
5135  
5140  
5145  
5150  
5155  
5160  
5165  
5170  
5175  
5180  
5185  
5190  
5195  
5200  
5205  
5210  
5215  
5220  
5225  
5230  
5235  
5240  
5245  
5250  
5255  
5260  
5265  
5270  
5275  
5280  
5285  
5290  
5295  
5300  
5305  
5310  
5315  
5320  
5325  
5330  
5335  
5340  
5345  
5350  
5355  
5360  
5365  
5370  
5375  
5380  
5385  
5390  
5395  
5400  
5405  
5410  
5415  
5420  
5425  
5430  
5435  
5440  
5445  
5450  
5455  
5460  
5465  
5470  
5475  
5480  
5485  
5490  
5495  
5500  
5505  
5510  
5515  
5520  
5525  
5530  
5535  
5540  
5545  
5550  
5555  
5560  
5565  
5570  
5575  
5580  
5585  
5590  
5595  
5600  
5605  
5610  
5615  
5620  
5625  
5630  
5635  
5640  
5645  
5650  
5655  
5660  
5665  
5670  
5675  
5680  
5685  
5690  
5695  
5700  
5705  
5710  
5715  
5720  
5725  
5730  
5735  
5740  
5745  
5750  
5755  
5760  
5765  
5770  
5775  
5780  
5785  
5790  
5795  
5800  
5805  
5810  
5815  
5820  
5825  
5830  
5835  
5840  
5845  
5850  
5855  
5860  
5865  
5870  
5875  
5880  
5885  
5890  
5895  
5900  
5905  
5910  
5915  
5920  
5925  
5930  
5935  
5940  
5945  
5950  
5955  
5960  
5965  
5970  
5975  
5980  
5985  
5990  
5995  
6000  
6005  
6010  
6015  
6020  
6025  
6030  
6035  
6040  
6045  
6050  
6055  
6060  
6065  
6070  
6075  
6080  
6085  
6090  
6095  
6100  
6105  
6110  
6115  
6120  
6125  
6130  
6135  
6140  
6145  
6150  
6155  
6160  
6165  
6170  
6175  
6180  
6185  
6190  
6195  
6200  
6205  
6210  
6215  
6220  
6225  
6230  
6235  
6240  
6245  
6250  
6255  
6260  
6265  
6270  
6275  
6280  
6285  
6290  
6295  
6300  
6305  
6310  
6315  
6320  
6325  
6330  
6335  
6340  
6345  
6350  
6355  
6360  
6365  
6370  
6375  
6380  
6385  
6390  
6395  
6400  
6405  
6410  
6415  
6420  
6425  
6430  
6435  
6440  
6445  
6450  
6455  
6460  
6465  
6470  
6475  
6480  
6485  
6490  
6495  
6500  
6505  
6510  
6515  
6520  
6525  
6530  
6535  
6540  
6545  
6550  
6555  
6560  
6565  
6570  
6575  
6580  
6585  
6590  
6595  
6600  
6605  
6610  
6615  
6620  
6625  
6630  
6635  
6640  
6645  
6650  
6655  
6660  
6665  
6670  
6675  
6680  
6685  
6690  
6695  
6700  
6705  
6710  
6715  
6720  
6725  
6730  
6735  
6740  
6745  
6750  
6755  
6760  
6765  
6770  
6775  
6780  
6785  
6790  
6795  
6800  
6805  
6810  
6815  
6820  
6825  
6830  
6835  
6840  
6845  
6850  
6855  
6860  
6865  
6870  
6875  
6880  
6885  
6890  
6895  
6900  
6905  
6910  
6915  
6920  
6925  
6930  
6935  
6940  
6945  
6950  
6955  
6960  
6965  
6970  
6975  
6980  
6985  
6990  
6995  
7000  
7005  
7010  
7015  
7020  
7025  
7030  
7035  
7040  
7045  
7050  
7055  
7060  
7065  
7070  
7075  
7080  
7085  
7090  
7095  
7100  
7105  
7110  
7115  
7120  
7125  
7130  
7135  
7140  
7145  
7150  
7155  
7160  
7165  
7170  
7175  
7180  
7185  
7190  
7195  
7200  
7205  
7210  
7215  
7220  
7225  
7230  
7235  
7240  
7245  
7250  
7255  
7260  
7265  
7270  
7275  
7280  
7285  
7290  
7295  
7300  
7305  
7310  
7315  
7320  
7325  
7330  
7335  
7340  
7345  
7350  
7355  
7360  
7365  
7370  
7375  
7380  
7385  
7390  
7395  
7400  
7405  
7410  
7415  
7420  
7425  
7430  
7435  
7440  
7445  
7450  
7455  
7460  
7465  
7470  
7475  
7480  
7485  
7490  
7495  
7500  
7505  
7510  
7515  
7520  
7525  
7530  
7535  
7540  
7545  
7550  
7555  
7560  
7565  
7570  
7575  
7580  
7585  
7590  
7595  
7600  
7605  
7610  
7615  
7620  
7625  
7630  
7635  
7640  
7645  
7650  
7655  
7660  
7665  
7670  
7675  
7680  
7685  
7690  
7695  
7700  
7705  
7710  
7715  
7720  
7725  
7730  
7735  
7740  
7745  
7750  
7755  
7760  
7765  
7770  
7775  
7780  
7785  
7790  
7795  
7800  
7805  
7810  
7815  
7820  
7825  
7830  
7835  
7840  
7845  
7850  
7855  
7860  
7865  
7870  
7875  
7880  
7885  
7890  
7895  
7900  
7905  
7910  
7915  
7920  
7925  
7930  
7935  
7940  
7945  
7950  
7955  
7960  
7965  
7970  
7975  
7980  
7985  
7990  
7995  
8000  
8005  
8010  
8015  
8020  
8025  
8030  
8035  
8040  
8045  
8050  
8055  
8060  
8065  
8070  
8075  
8080  
8085  
8090  
8095  
8100  
8105  
8110  
8115  
8120  
8125  
8130  
8135  
8140  
8145  
8150  
8155  
8160  
8165  
8170  
8175  
8180  
8185  
8190  
8195  
8200  
8205  
8210  
8215  
8220  
8225  
8230  
8235  
8240  
8245  
8250  
8255  
8260  
8265  
8270  
8275  
8280  
8285  
8290  
8295  
8300  
8305  
8310  
8315  
8320  
8325  
8330  
8335  
8340  
8345  
8350  
8355  
8360  
8365  
8370  
8375  
8380  
8385  
8390  
8395  
8400  
8405  
8410  
8415  
8420  
8425  
8430  
8435  
8440  
8445  
8450  
8455  
8460  
8465  
8470  
8475  
8480  
8485  
8490  
8495  
8500  
8505  
8510  
8515  
8520  
8525  
8530  
8535  
8540  
8545  
8550  
8555  
8560  
8565  
8570  
8575  
8580  
8585  
8590  
8595  
8600  
8605  
8610  
8615  
8620  
8625  
8630  
8635  
8640  
8645  
8650  
8655  
8660  
8665  
8670  
8675  
8680  
8685  
8690  
8695  
8700  
8705  
8710  
8715  
8720  
8725  
8730  
8735  
8740  
8745  
8750  
8755  
8760  
8765  
8770  
8775  
8780  
8785  
8790  
8795  
8800  
8805  
8810  
8815  
8820  
8825  
8830  
8835  
8840  
8845  
8850  
8855  
8860  
8865  
8870  
8875  
8880  
8885  
8890  
8895  
8900  
8905  
8910  
8915  
8920  
8925  
8930  
8935  
8940  
8945  
8950  
8955  
8960  
8965  
8970  
8975  
8980  
8985  
8990  
8995  
9000  
9005  
9010  
9015  
9020  
9025  
9030  
9035  
9040  
9045  
9050  
9055  
9060  
9065  
9070  
9075  
9080  
9085  
9090  
9095  
9100  
9105  
9110  
9115  
9120  
9125  
9130  
9135  
9140  
9145  
9150  
9155  
9160  
9165  
9170  
9175  
9180  
9185  
9190  
9195  
9200  
9205  
9210  
9215  
9220  
9225  
9230  
9235  
9240  
9245  
9250  
9255  
9260  
9265  
9270  
9275  
9280  
9285  
9290  
9295  
9300  
9305  
9310  
9315  
9320  
9325  
9330  
9335  
9340  
9345  
9350  
9355  
9360  
9365  
9370  
9375  
9380  
9385  
9390  
9395  
9400  
9405  
9410  
9415  
9420  
9425  
9430  
9435  
9440  
9445  
9450  
9455  
9460  
9465  
9470  
9475  
9480  
9485  
9490  
9495  
9500  
9505  
9510  
9515  
9520  
9525  
9530  
9535  
9540  
9545  
9550  
9555  
9560  
9565  
9570  
9575  
9580  
9585  
9590  
9595  
9600  
9605  
9610  
9615  
9620  
9625  
9630  
9635  
9640  
9645  
9650  
9655  
9660  
9665  
9670  
9675  
9680  
9685  
9690  
9695  
9700  
9705  
9710  
9715  
9720  
9725  
9730  
9735  
9740  
9745  
9750  
9755  
9760  
9765  
9770  
9775  
9780  
9785  
9790  
9795  
9800  
9805  
9810  
9815  
9820  
9825  
9830  
9835  
9840  
9845  
9850  
9855  
9860  
9865  
9870  
9875  
9880  
9885  
9890  
9895  
9900  
9905  
9910  
9915  
9920  
9925  
9930  
9935  
9940  
9945  
9950  
9955  
9960  
9965  
9970  
9975  
9980  
9985  
9990  
9995  
10000  
10005  
10010  
10015  
10020  
10025  
10030  
10035  
10040  
10045  
10050  
10055  
10060  
10065  
10070  
10075  
10080  
10085  
10090  
10095  
10100  
10105  
10110  
10115  
10120  
10125  
10130  
10135  
10140  
10145  
10150  
10155  
10160  
10165  
10170  
10175  
10180  
10185  
10190  
10195  
10200  
10205  
10210  
10215  
10220  
10225  
10230  
10235  
10240  
10245  
10250  
10255  
10260  
10265  
1



polymer. This method includes not only the charge density contributions from the copolymerization reactants, but also includes the anionicity that results from hydrolysis of the amide group that occurs during the polymerization reaction. The hydrolysis reaction can convert a significant number of the amide groups to carboxyl groups, as much as 7% for low charge density polymers (as in the case of the Floerger AN 905 PWG polymer). In our experience, it is the total anionicity that dictates how well a given polymer performs in this process application. "Polymer charge density" or "polymer anionicity" herein is defined as the total charge density, which includes charge resulting from the copolymerization reaction (i.e. theoretical charge density), plus the charge contributions originating from hydrolysis.

The anionic features of the polymeric flocculants can be imparted by any appropriate chemical constituents present on the polymer. Anionic polymers containing carboxyl, carboxymethyl, phosphate and sulfate functionalities, for example, are particularly well suited for the invention. Examples of preferred flocculants according to this invention include anionic polyacrylamides and anionic polysaccharides.

Polyacrylamides with carboxyl groups represent one preferred class of flocculants for use in this invention. Those having polymer charge densities between about 1 and 30%, more preferably between about 1 and 20%, and most preferably between about 5 and 12% at neutral pH, have been found most suitable.

A flocculant used in this invention can have essentially any molecular weight provided it does not adversely effect the desired flocculation of the insoluble contaminants. Higher molecular weights may provide improved flocculation performance. Typically, the average molecular weight will be in the range from about 100,000 to 50,000,000 or more. Preferably it is greater than about 100,000, more preferably greater than about 1,000,000, and most preferably greater than about 10,000,000, provided this does not adversely effect the desired flocculation reaction.

In a preferred embodiment of this invention, a ST protein suspension is clarified by the above approach. U.S. Patent No. 5,182,369, the disclosure of which is incorporated herein by reference, describes the selective precipitation of somatotropin oligomers and other contaminants

as a means for the recovery of highly purified somatotropin monomer from a refold solution. This is accomplished by reducing the pH of the refold solution from the high level employed for the naturation step (usually in excess of pH 10) down to an optimum pH end point value that is generally in the range of about 4.0 to about 6.5.

5           The particles that result from the precipitation step are generally small and do not readily settle. Under gravitational conditions, a significant quantity of the solids stay suspended for an indefinite period of time such that a liquid/solid interface is not visually observable. The characteristically slow sedimentation velocity of this suspension makes clarification using centrifugation difficult. Depth filtration has also been used to clarify this process stream, but  
10 there are also disadvantages associated with this type of separation, such as the expense of additional raw materials and waste handling.

Although additional purification steps may be desired, for example by chromatographic or other methods, the ST monomer recovered in accordance with this invention is substantially free of residual bacterial protein and other contaminants and is suitable for administration to the target animal by injection or implantation without further purification. The purified product preferably contains less than about 10%, more preferably less than about 5% of somatotropin oligomers. Such oligomers are biologically inactive and do not result in any adverse reaction when administered to the target animal.

15  
20           The mixing conditions used in the flocculation process can effect several parameters, as is known in the art. These parameters include the sedimentation rate of the flocs and the clarity of the supernatant. The preferred parameters with respect to mixing often represent a compromise designed to provide the optimal flocculation conditions possible in the face of frequently conflicting requirements of the individual parameters. For example, and increase in supernatant clarity is typically accompanied by a decrease in the sedimentation rate of  
25 the flocs.

Mixing of the polymer solution into the process pool can depend on not only the size, type, location, orientation and rotational speed of the mixer, but also can depend on the rate at which the polymer solution is delivered. So for a given impeller size, two parameters that can be

manipulated include mixer rotational speed and the rate of polymer solution addition. The skilled individual in the art can readily determine the optimal conditions in this regard for a given product and process configuration.

5 A balance must be reached between agitation intensity and shear rate. The design of the mixing system used in accordance with this invention can assist in achieving this balance. A mixing impeller design that has a high pumping capacity and low shear rate is well suited for this application. For example, the Chemineer HE-3 high efficiency impeller (Chemineer - Dayton, OH) or the Lightnin A-310 high efficiency mixer impeller (Lightnin Mixers - Rochester, NY) are suitable choices for this type of operation.

10 The mixer should also be capable of keeping the flocs suspended during the operation. The minimum level of agitation necessary in a flocculation system is bound by the fluid velocity needed to keep the formed flocs in a "just suspended state", so that solids are not allowed to settle at the bottom of the tank (Gates et al., 1976).

15 As alluded to above, two desirable features of the flocculation system described herein are good mixing and low shear. Tank design can to some extent address both of these concerns. Low shear can be met by minimizing the number of surfaces that cause shear by either creation of fluid motion or by surface particle interaction. Proper tank geometry and flocculant feed location will help ensure uniform mixing of the flocculant into the process fluids.

20 The following examples are included to demonstrate preferred embodiments of the invention. It should be appreciated by those of skill in the art that the techniques disclosed in the examples which follow represent techniques discovered by the inventor to function well in the practice of the invention, and thus can be considered to constitute preferred modes for its practice. However, those of skill in the art should, in light of the present disclosure, appreciate that many changes can be made in the specific embodiments which are disclosed and still obtain  
25 a like or similar result without departing from the spirit and scope of the invention.

## EXAMPLE 1

The objective of the experiments described herein was to survey and identify flocculating agents that would aggregate the suspended particles from an aqueous suspension containing soluble and insoluble components. In this example, aqueous solutions of the flocculating agent were made up and added to the ST protein suspension. The flocculating agents were blended into this precipitate suspension using an overhead mixer as described in the Lab-Scale Flocculation Procedure # 1. Upon completion of the flocculating agent addition, the mixer was stopped and the flocculated solids allowed to settle. The sedimentation velocity and supernatant clarity were then measured as described below:

### Lab-Scale Flocculation Procedure #1

The experimental lab-scale flocculation equipment was assembled using a 4-liter glass beaker, a Lightnin® LabMaster™ mixer with a 6.35 cm diameter A-310 high efficiency mixer impeller, an Ismatec peristaltic pump (model mv-ge) and a container for the polymer solution. The 4-liter glass beaker was baffled to enhance vertical mixing in the system. The 4-liter glass beaker was designed with glass indentions that ran from 1 cm above the bottom of the beaker to just below the lip of the beaker. The pump tubing size was selected based on the required feed flow rates for the polymer solution and a 1/8" stainless steel dip tube was used to feed the polymer solution. The outlet of this tube was located just below and just outside the impeller perimeter so as to ensure efficient mixing of the polymer into the process pool. The Ismatec pump was calibrated to obtain the desired polymer feed flow rate. A strip of tape with 1 cm increments marked on it was applied in a vertical position on the side of the 4-liter beaker to provide a means of measuring the level of the solid-liquid interface. The impeller was positioned approximately 3 to 4 cm above the bottom of the 4-liter beaker. The impeller position was fixed throughout a set of experiments by aligning the mixer chuck with a mark on the impeller shaft.

Two liters of the ST protein suspension were transferred to the 4-liter beaker, and then the mixer was turned on with the rotational speed set at 169 rpm. One liter of the polymer solution was pumped into the process pool at a constant rate over a 30 minute time frame. After delivery

of the polymer solution into the process pool, the mixer was allowed to mix the contents for an additional 60 seconds and then the mixer was turned off. The settling rate of the flocculated particles and the supernatant absorbance were then measured. The results for a variety of flocculating agents are shown in Table 1.

## 5 Measurements:

10 Sedimentation Velocity: Upon completion of the polymer solution addition to the precipitate solution, the mixer was stopped, and the flocculated solids were allowed to settle. After the initial momentum of the fluids diminished, the flocs began to settle and a visibly discrete solid-liquid interface is formed. Sedimentation velocity was measured by timing the change in the solid-liquid interface level. The timing started immediately after the mixer was stopped and the timer was stopped at a designated solid-liquid interface level in the glass vessel. The velocity at the interface represents the velocity of the slowest moving particles, which have sufficient concentration to be visible. The slowest moving particles are rate limiting with respect to cycle time and therefore are the most interesting in terms of characterizing the flocculation process. The accuracy of this measurement technique is dependent on the formation of a discernible solid-liquid interface and on the sedimentation velocity. The formation of large particles creates a situation where sedimentation rates are fast and where a discernible solid-liquid interface is not formed until the solids have settled to the bottom of the container. In such cases the measurement of sedimentation velocity is difficult and will likely have a larger margin of error.

20 Absorbance @ 700 nm: Measurement of the absorbance at 700 nm ( $A_{700}$ ) is a means of quantitatively assess the supernatant clarity achieved from a given treatment or process. The absorbance of a 1 ml sample of the supernatant after treatment was measured using a Hewlett Packard UV-VIS spectrophotometer (model 8453) and a cuvette with a 1 cm path length. The absorbance value obtained using this method can be converted to a transmittance value as per equation 14:

$$\%T = 100 \times \exp(-A_{700nm})$$

Intuitively, transmittance is easy to relate to clarity of the supernatant since it is directly related to clarity, whereas absorbance is inversely related to clarity. The clarification process of this invention results in a solution having  $A_{700}$  less than 0.5, preferably less than 0.1.

**Table 1: Survey of flocculating agents**

Flocculant Trade Name	Manufacturer	Characteristics					Effectiveness (+/-)
		Composition	MW <sub>ave</sub> (Daltons)	Anion/ Cationic	Charged Group	Theoretical Charge Density	
STA-COTE H-44	A. E. Staley Mfg. Co.	Corn Starch					-
STA-JEL 140	A. E. Staley Mfg. Co.	Potato Starch	1.5 M	Anionic	phosphate		+
Floerger FA 920	Chemtall, Inc.	Polyacrylamide	15-20 M	Nonionic		0%	-
Floerger AH 912	Chemtall, Inc.	Polyacrylamide	15-20 M	Anionic	hydrolyzed homopolymer	~2%	+
Floerger AN 903 PWG/GR	Chemtall, Inc.	Polyacrylamide	15-20 M	Anionic	carboxyl	1.5%	++
Floerger AN 905 PWG/GR/VHM	Chemtall, Inc.	Polyacrylamide	15-20 M	Anionic	carboxyl	3%	++
Floerger AN 910 PWG	Chemtall, Inc.	Polyacrylamide	15-20 M	Anionic	carboxyl	10%	++
Floerger AN 923 PWG	Chemtall, Inc.	Polyacrylamide	15-20 M	Anionic	carboxyl	20%	+
Floerger AN 934 PWG	Chemtall, Inc.	Polyacrylamide	15-20 M	Anionic	carboxyl	30%	+
Agefloc WT601	Ciba	Polyacrylamide	15-20 M	Anionic	carboxyl	10%	+
Agefloc WT611	Ciba	Polyacrylamide	15-20 M	Anionic	carboxyl	10%	++
CYSEP 2411	CYTEC, Inc.		~ 6 M	Cationic			-
CYSEP 2709	CYTEC, Inc.		~ 4 M	Cationic			-
CYSEP 329	CYTEC, Inc.	Polyamine	~ 50 K	Cationic			-
CYSEP 349	CYTEC, Inc.	Polyamine	~ 250 K	Cationic			-
CYSEP 4022	CYTEC, Inc.	PolyDADM	~ 50 K	Cationic			-
CYSEP 615	CYTEC, Inc.	Polyamine	~ 30 K	Cationic			-
Methocel A4M PREM	Dow Chemical Co.	Methcellulose		Nonionic			-
Methocel E4M PREM	Dow Chemical Co.	Hydroxypropyl Methocellulose		Nonionic			-
Methocel F4M PREM	Dow Chemical Co.	Hydroxypropyl Methocellulose		Nonionic			-
Methocel K100M PREM	Dow Chemical Co.	Hydroxypropyl Methocellulose		Nonionic			-
Chitosan	SIGMA	Chitosan	~ 370 K				-
Prestol 2500 TR	Stockhausen, Inc.	Polyacrylamide	15-20 M	Nonionic		0%	-
Prestol 2505	Stockhausen, Inc.	Polyacrylamide	15-20 M	Anionic	Carboxyl	5%	+
Prestol 2515 TR	Stockhausen, Inc.	Polyacrylamide	15-20 M	Anionic	Carboxyl	15%	++
Prestol 2520	Stockhausen, Inc.	Polyacrylamide	15-20 M	Anionic	Carboxyl	20%	+

"+" = clear supernatant (%Transmittance  $\geq$  70%), slow sedimentation (velocity < 500 cm/hr)

"++" = clear supernatant (%Transmittance  $\geq$  95%), fast sedimentation (velocity  $\geq$  500 cm/hr); "NT" = not tested;

"-" = no visible formation of a solid-liquid interface; "M" = millions; "K" = thousands; MW = molecular weight

These results show that only anionic flocculating agents produced detectable signs of flocculation as indicated by the sedimentation velocity of the flocculated solids and by the clarity (reported as % transmittance) of the aqueous phase (supernatant). Flocculating agents producing

flocs that exhibited sedimentation velocities greater than 500 cm/hr and supernatants with % transmittance greater than 95% were given a rating of “++”. Flocculating agents producing flocs that exhibited sedimentation velocities less than 500 cm/hr and supernatants with % transmittance greater than 70% were given a rating of “+”. Flocculating agents producing no visible signs of flocculation, as indicated by the lack of formation of a liquid-solid were given a rating of “-”. It should be noted that the reported charge densities for the respective flocculating agents are theoretical charge densities, and do not represent the total charge density.

## EXAMPLE 2

A range of flocculating agent concentrations was evaluated using Lab-Scale Flocculation Procedure #2. A high molecular weight (roughly 16 million daltons) acrylamide polymer with anionic character was used in this evaluation. The three primary measurements of interest in these experiments were supernatant absorbance, total protein concentration and bST concentration. Table 2 displays the data obtained using this procedure. The final polymer concentration referred to herein represent the total polymer concentration after all of the polymer solution had been added to the ST protein suspension.

### Lab-Scale Flocculation Procedure #2

1. Make up the flocculant solution at room temperature to the desired concentration (i.e. 200 ppm) 1 hour in advance of the flocculation experiments.
2. Label 3 polycarbonate conical-bottom test tubes (w/ screw caps). Each polymer concentration was evaluated in triplicate.
3. Pipette the ST protein suspension into each test tube. Pipette the corresponding volume of the polymer solution into the ST protein suspension and then immediately cap the test tube and invert 20 times.
4. Place the test tubes into a bucket type centrifuge and spin at room temperature for 5 minutes at 1000 rpm (Centrifuge type: Sorvall RC5C, Rotor size: SS34, RCF ~ 119).

5. After centrifugation, gently remove the cap and using a transfer pipette pull the top 6 ml (of 10 ml total) of supernatant. Draw liquid near the liquid/air interface while trying to avoid floating solids. Transfer the supernatant to a labeled test tube.

6. Measure the absorbance of the supernatant at 700 nm. Then dilute 1 ml of the supernatant in 5 ml of 1% acetic acid. Measure the absorbance at 280 nm and use the supernatant for HPLC analysis to determine bST concentration.

Repeat steps 2 through 6 for a range of polymer concentrations.

#### Measurements:

Total Protein Concentration: The total protein concentration was calculated from a measured absorbance at a wavelength of 278 nm on a diluted supernatant sample. Solutions of bST have an extinction coefficient (extinction coefficient,  $\epsilon$ , or molar absorptivity is defined as the constant of proportionality in Beer's law:  $A = \epsilon bc$ , where  $A$  is absorbance,  $b$  is pathlength, and  $c$  is the molarity of the absorbing species) of 0.651. The supernatant sample was diluted in 1% acetic acid to obtain a total protein concentration near 1.5 grams/liter. The use of acetic acid as the dilution buffer serves to dissolve any suspended solids not removed in the treatment. (The suspended solids not removed in the treatment are considered part of the supernatant. These solids must be dissolved so they do not interfere with the absorbance measurement at 278 nm.) For this measurement, a 1% acetic acid solution was used as the zero reference.

bST Concentration: A portion of the diluted sample from the total protein concentration measurement was placed into a vial and analyzed using reverse phase HPLC.

The first row of data in Table 2 are the data associated with a sample of the ST protein solution that was filtered through a 0.2 micron syringe filter. The syringe filter sample serves as a means to assess the interaction, if any, of the flocculating agent with the soluble ST protein. The second row of data in Table 2 are data associated with the "control". The control is a sample of ST protein suspension to which flocculating agent was not added.



**Table 2: Effect of Flocculant Polymer Dosage on the Centrifugal Clarification of pH Adjusted ST Protein Suspension**

Polymer Concentration (ppm)	Supernatant Absorbance (A700)	Supernatant Protein Concentration (g/L)*	bST Concentration (g/L)*	Yield (%)
0	0.01	8.50	<b>7.72</b>	100%
0	1.67	13.48	8.71	113%
1	1.65	12.53	8.62	112%
10	1.32	8.97	8.12	105%
50	0.15	8.13	8.04	104%
100	0.06	7.87	8.15	106%
500	0.30	6.93	6.75	87%
1000	0.67	2.35	3.18	41%

\*Note: Protein and bST concentrations are corrected for dilutions (i.e. they reflect predilution values).

The data show that the total protein concentration decreased, as did the supernatant absorbance, upon addition of the flocculating agent for final flocculating agent concentrations greater than 1 ppm, when compared to the control, while the soluble ST concentration remained relatively constant for samples that were treated with final flocculating agent concentrations less than 500 ppm. From these data we can conclude the following: 1) The flocculating agent did not exhibit significant interaction with the soluble ST protein so as to remove the soluble protein from solution for final flocculating agent concentrations less than 500 ppm. 2) The decrease in supernatant absorbance and the decrease in the total protein concentration indicates that the addition of the flocculating agent to the solution containing soluble ST protein and precipitated contaminants improved the sedimentation properties of the solids in this system. The increase in sedimentation performance is directly related to the size of the particles in this system, which can be attributed the flocculation activity of the flocculating agent. Flocculation as measured by monitoring the total protein concentration was observed to occur for flocculating agent concentrations as low as 1 ppm. Flocculating agent concentrations above 50 ppm yielded purities at or above 95%, but the bST concentration began to decrease for flocculating agent concentrations greater than about 100 ppm by this procedure. Thus, although a preferred range of polymer concentrations based on these experiments was between about 5 and 100 ppm, the invention is operable at up to 1000 ppm or more of polymer. The observed decrease in bST concentration, at polymer concentrations greater than or equal to 500 ppm, may be an indication that the flocculant agent was interacting with bST and pulling it down with the other solids.

### EXAMPLE 3

Polyacrylamide flocculating agents manufactured by Chemtall Inc. (Riceboro, GA) were evaluated. The trade name for the Chemtall flocculating agent tested in this example is Floerger AN 905 GR. Six lab-scale flocculation experiments were run each using 2 liters of ST protein suspension and 1 liter of flocculant solution. The flocculant solution was added to a mixed ST protein suspension at a constant flow rate over an approximate time frame of 30 to 40 minutes. The mixing conditions were held constant for this set of experiments so that the effect of the polymer concentration on the flocculation performance could be assessed. The sedimentation velocity was approximated by measuring the time it took for the liquid/solid interface to reach a fixed level in the glass beaker after the over-head mixer was turned off. Supernatant absorbance was measured using a spectrophotometer. Table 3 illustrates flocculation performance data for a range of Floerger An 905 PWG flocculant concentrations.

Table 3 shows the results generated when polymeric acrylamide flocculant with an average molecular weight of about 16 million daltons and a charge density of 10% (Floerger AN 905 GR [Lot # 8S36839]) was added to the ST protein suspension to obtain final polymer concentrations ranging from 5 to 35 ppm. The mixer speed was 169 rpm, the tank diameter was 15.9 cm, temperature was about 20 °C, mixer impeller pump capacity was 24 liters/minute, impeller diameter was 6.35 cm, and impeller tip speed was 56 cm/sec, with center mixing in a baffled container.

The data show that the sedimentation velocity increased with an increase in the final polymer concentration, while the supernatant transmittance decreased with an increase in final polymer concentration. The increase in sedimentation velocity is an expected outcome in that more polymer is available to form larger aggregates. The larger aggregates will settle faster as a result of their larger size.

**Table 3: Effect Polymer Concentration on the Flocculation Performance**

Exp. #	bST Volume (L)	[AN905] <sub>0</sub> (ppm)	Polymer Volume Added (L)	[AN905] <sub>f</sub> (ppm)	Polymer Addition Time (min)	Flocculation Performance		
						Initial Sedimentation Velocity (cm/hr)	Supernatant Absorbance ( $\lambda = 700$ nm)	Supernatant Transmittance (%)
1	2	15	1	5	30	77	0.011	99%
2					30	78	0.016	98%
3					30	85	0.003	100%
					Average =	80		99%
					Std. Dev. =	4		0.63%
4	2	45	1	15	30	989	0.026	97%
5					30	1147	0.022	98%
6					30	1042	0.022	98%
					Average =	1059		98%
					Std. Dev. =	80		0.25%
7	2	75	1	25	30	1205	0.040	96%
8					31	1266	0.039	96%
9					30	1456	0.031	97%
					Average =	1309		96%
					Std. Dev. =	131		0.45%
10	2	105	1	35	30	1603	0.047	95%
11					30	1767	0.052	95%
12					30	2231	0.038	96%
					Average =	1685		96%
					Std. Dev. =	116		0.70%

Sedimentation velocity data point was not included in the average or std. dev. because the solids did not suspend at during the post mixing stage of the run. The flocculated contents were mixed for 1 minute after completion of the polymer addition.

#### EXAMPLE 4

Data on the effect of the polymer charge density on the flocculation performance

is listed in Table 4. Two final polymer concentrations (15 ppm and 30 ppm) were evaluated

- 5 using Chemtall flocculant polymers. The polymers used in this set of experiments had theoretical charge densities ranging from 0% to 20%. The reaction conditions were as follows: ST protein suspension volume = 2 liters, mixer size, capacity, and speed as in EXAMPLE 3, tank diameter = 15.9 cm, 1/8" dip tube, volume of polymer added = 1 liter resulting in a 50% increase in volume.

**Table 4: Summary of Chemtall Flocculant Evaluations**

Exp. #	Flocculant	Theoretical Anionicity (%)	Measured Anionicity (%)	[F] <sub>0</sub> (ppm)	[F] <sub>r</sub> (ppm)	Actual polymer addition time (min)	Pump speed setting	Flocculation Performance		
								Initial sedimentation velocity (cm/hr)	Supernatant absorbance ( $\lambda = 700$ nm)	Transmittance (%)
1	Floerger FA 920	0	1	45	15	31	500	0	1.297	27%
2	Floerger AN 903 PWG	1.5	6	45	15	32	500	785	0.053	95%
3	Floerger AN 905 PWG	3	10	45	15	32	510	947	0.034	97%
4	Floerger AN 910 PWG	10	14	45	15	32	520	275	0.017	98%
5	Floerger 923 AN PWG	20	27	45	15	31	530	23	0.026	97%
6	Floerger FA 920	0	1	90	30			0	1.313	27%
7	Floerger AN 903 PWG	1.5	6	90	30	31	540	1284	0.050	95%
8	Floerger AN 905 PWG	3	10	90	30	31	540	1657	0.024	98%
9	Floerger AN 910 PWG	10	14	90	30	32	540	391	0.002	100%
10	Floerger 923 AN PWG	20	27	90	30		550	25	0.026	97%

The data in Table 4 was generated using Chemtall flocculating agents. These data were obtained using the Lab-Scale Flocculation Procedure #1. The sedimentation velocity and the supernatant transmittance data were plotted against the total polymer anionicity (as measured by the manufacturer), see Figure 1. Table 4 lists data of the flocculation performance parameters with the actual measured polymer anionicity. Only anionic polymers manufactured using the acrylate and acrylamide polymerization were used to produce the data in Table 4 (as opposed to using anionic polymers that were made by hydrolysis). The theoretical anionicity is calculated by the amount of acrylate used in the copolymerization reaction. The total anionicity includes the theoretical anionicity, plus the anionicity resulting from hydrolysis.

The data in Table 4 show that polymeric flocculants that have an anionicity of 1% or lower did not exhibit detectable signs of flocculation as evidenced by the lack of formation of a visible solid-liquid interface at static conditions. The data points corresponding to a 1% anionicity represent the performance of a polymer that was manufactured as a nonionic polymer, that is, a polymer without carboxyl groups. The 1% anionic character of this polymer is the

result of acid hydrolysis that occurred during the manufacturing process. The anionic character of a polymer containing carboxyl groups is significantly reduced at a low pH. This fact implies that a 1% anionic charge on the polymer at neutral pH (pH ~ 7.0) is essentially neutralized (i.e. charge density  $\approx 0$ ) at pH 4.5. For purposes of this discussion, the polymer associated with the 1% anionicity data points is considered as having a nonionic character. Formation of a solid-liquid interface did not occur at static conditions which suggests this polymer did not significantly affect the size of the precipitate particles in the suspension. Measurement of the solution clarity was essentially measurement of the diluted suspension since a separation of solids from the liquid did not occur, which explains the low transmittance value indicated in Table 4.

Flocculations performed using polymers with charge densities greater than 1% exhibited a discernible solid-liquid interface. The data clearly show that there exists an optimum anionicity for flocculation of the solids in this system. The 10% anionic polymer exhibited the fastest sedimentation velocity, which indicates that the optimum anionicity is near 10%. Sedimentation velocities obtained using polymers with charge densities less than or greater than the 10% anionic polymer were measurably slower than the velocity obtained when using the 10% anionic polymer, while the supernatant clarity appears to be relatively consistent for all flocculations performed using polymers with anionic character.

A polymeric flocculant with a charge density less than the optimal charge density (i.e. the 6% anionic polymer) exhibited a slower sedimentation velocity than the velocity achieved when using the 10% anionic polymer. The data indicate that charge densities less the optimal charge density result in less efficient flocculation of the solids.

## EXAMPLE 5

The effect of the pH of the pH-adjusted ST protein suspension on flocculation performance was evaluated. The total protein concentration of refold concentrate was adjusted to approximately 20 g/L using purified water. The pH of the diluted refold solution was adjusted

using 5% acetic acid to 4 different pH endpoints. At each pH a 2 liter sample of the resulting suspension was taken and further mixed for at a minimum of 1 hour using a stir bar and stir plate. Each of the suspension samples were flocculated using the Lab-Scale Flocculation Procedure #1 with a 1 liter aliquot of a 75 ppm solution of Floerger AN 905 PWG flocculating agent (Chemtall Inc., 10% anionicity). The dip tube diameter was 1/8", the mixer rotational speed was 169 rpm. The pH adjustment and the flocculations were run at room temperature (~ 23° C). The pH data and the flocculation performance data is listed in the Tables 5 and 6 below. These data indicate that for the Floerger AN 905 PWG polyacrylamide flocculant (F905), the preferred pH of the pH-adjusted ST protein suspension is between about 4 and 5, more preferably between about 4.2 and 4.6.

**Table 5: Flocculation Dependence on pH**

Exp. #	pH	Pump speed setting	Q <sub>P</sub> (ml/min)	Polymer addition time (min)	Flocculation Performance		
					Initial Sedimentation velocity (cm/hr)	Supernatant Absorbance (700 nm)	%T
1	4.73	600	37	27	504	0.032	97
2	4.60	600	37	27	932	0.053	95
3	4.50	600	37	27	1219	0.026	97
4	4.40	600	37	27	1477	0.023	98
5	4.30	600	37	27	1567	0.025	98

**Table 6: Flocculation Dependence on pH**

Exp. #	pH	Pump speed setting	Q <sub>P</sub> (ml/min)	Polymer addition time (min)	Flocculation Performance		
					Initial Sedimentation velocity (cm/hr)	Supernatant Absorbance (700 nm)	%T
1	4.80	810	31	32	104	0.014	99
2	4.60	810	31	32	686	0.016	98
3	4.50	810	31	32	960	0.013	99
4	4.40	810	31	32	1527	0.016	98
5	4.30	810	31	32	1292	0.009	99

The data in Tables 5 and 6 show that flocculation performance, as measured by sedimentation velocity, is strongly dependent on solution pH. These data are represented in graphical form in Figures 2 and 3. The sedimentation velocity significantly increases as the pH decreases from 4.8 to 4.3. The clarity, as measured by % transmittance, of the supernatant in this pH range is relatively constant. The data in Tables 5 and 6 also indicate that for high molecular

weight polymer of about 16,000,000 daltons, with a measured charge density of about 10%, the flocculation performance as measured by the sedimentation velocity is optimal near about pH 4.40.

5

## EXAMPLE 6

A modified version of Lab-Scale Flocculation Procedure #2 was used to evaluate starch and cellulose polymers as flocculating agents. Essentially the same procedure was used as above except for the process volumes, beaker sizes and polymer addition times. In this experiment 1 liter of polymer was added to 1 liter of ST protein suspension over a 1 to 2 hour time frame.

10

Starch polymers were provided by A. E. Staley Manufacturing Company (Decatur, IL). One potato starch (Sta Jel 140, prejelatenized potato starch) were evaluated. The potato starch exhibited visual signs of particle aggregation under the conditions tested.

009399215  
10/28/01 10:27:01

20

Sta Jel 140 is a prejelatenized potato starch, a starch which is primarily composed of amylose and amylopectin polymers. Amylose is linear polymer, whereas the amylopectin is a highly branched polymer. Most starches contain 18 to 28% amylose, with potato starch exhibiting amylose content at the lower end of this range (Wurzburg). Potato starch contains approximately 20% amylose and 80% amylopectin, and has an approximate average molecular weight of 1.5 million daltons. This ratio of amylose to amylopectin suggests that potato starch is a highly branched polymer. Unmodified potato starch naturally contains 0.08% phosphorous, which gives this starch its anionic characteristics.

25

Potato feed starch concentrations of 0.02% (wt/wt) and 0.2% (wt/wt) were tested to obtain final starch polymer concentrations of 0.01% (wt/wt) and 0.1% (wt/wt), respectively. In both cases the Sta Jel 140 potato starch exhibited visual signs of flocculation. The solids settled at gravitational conditions leaving a nearly crystal clear supernatant. The formed flocs were smaller than the flocs formed using the Floerger AN 905 PWG polymer, which resulted in slower sedimentation of the solids. Measurement of the total soluble protein indicated that there was no detectable loss in soluble protein concentration upon treatment of the ST protein suspension with

the potato starch. This is evidence that the potato starch did not interact with the soluble bST protein.

All of the compositions and methods disclosed and claimed herein can be made and executed without undue experimentation in light of the present disclosure. While the compositions and methods of this invention have been described in terms of preferred embodiments, it will be apparent to those of skill in the art that variations may be applied in the steps or in the sequence of steps of the method described herein without departing from the concept, spirit and scope of the invention. More specifically, it will be apparent that certain agents which are both chemically and physiologically related may be substituted for the agents described herein while the same or similar results would be achieved. All such similar substitutes and modifications apparent to those skilled in the art are deemed to be within the spirit, scope and concept of the invention as defined by the appended claims.